



I.V. POLE STANDS AND TOOL-FREE ACCESSORIES



Medline has an I.V. Pole Stand to meet all your needs and budgets — from standard duty to fully-featured heavy-duty models.

PREMIER HEAVY-DUTY I.V. POLE

WITH QUICK RELEASE CASTERS


Ideal for triple-channel pump or any application that would require an extra stable base

- Patent-pending 3" (7.6 cm) premium Quick Release Casters help keep pole in service longer
- Heavy-duty 24" (61 cm) diameter low center of gravity hammer-tone cast iron base with 6 wheels help resist tipping, weighs 19.7 lbs. (9 kg)
- Rubber bumpers at end of legs help to protect walls and other objects from scuffing
- D-tube insert helps ensure stability of pole without the fear of base spinning
- Stainless steel telescoping pole adjusts from 48" (122 cm) compressed to 82" (208 cm) fully extended
- Four durable ram hooks

TOOL-FREE ACCESSORIES FOR MDS80600 MODEL

URINE BAG HOOK

Item#: MDS80600UBH*
Pkg: 10/cs



STAND ASSIST WITH TRAY

Item#: MDS80600SAT*
Pkg: 6/cs



GRAB BAR

Item#: MDS80600GB*
Pkg: 6/cs



BASKET

Item#: MDS80600B*
Pkg: 2/cs



OXYGEN TANK HOLDER

Item#: MDS80600OTH*
Pkg: 2/cs



CASTERS

Item#: MDS80600C*
Pkg: 2/bx



ADD VERSATILITY

* Patent Pending

ONE OF THE MOST STABLE BASES AVAILABLE!



pump is not included

MDS80600SAT (optional)



MDS80600B (optional)

MDS80600UBH (optional)

ITEM# MDS80600
PKG. 2/CS



TO ORDER: 1-800-MEDLINE


HEAVY-DUTY ALUMINUM

I.V. POLE STAND

Ideal for pumps, multiple bags, or any application which requires an extra stable base

- Constructed of corrosion-resistant aluminum with heavy-duty 18 lb. (8.2 kg) cast iron five-legged base with low center of gravity to help ensure stability
- Base diameter 21.5" (54.6 cm)
- Accommodates two pumps up to 14 lbs. (6.3 kg) each and two 3-liter bags of fluid for versatility
- Friction grip lock makes adjusting easy
- Height adjustable from 46.75" to 81" (119 cm to 205 cm)
- Four durable ram hooks hold items securely

ITEM# MDS80500
PKG. 2/CS




DELUXE

I.V. POLE STAND

Standard duty model featuring an extra-high height range

- Constructed of corrosion-resistant aluminum with five-legged, 23.5" diameter (59.7 cm) cast aluminum base weighing 10 lbs. (4.5 kg)
- 2" (5 cm) heavy-duty casters provide smooth, quiet movement
- Accommodates two pumps up to 14 lbs. (6.3 kg) and two 3-liter bags of fluid for versatility
- Microphone-style knob makes adjustments easy
- Height adjustable from 56.75" to 96" (144 cm to 244 cm)
- Four ram hooks hold items securely

ITEM# MDS80494
PKG. 2/CS



STANDARD

I.V. POLE STAND

Value-positioned to meet basic needs

- Large 23" (58 cm) diameter, four-legged base provides excellent stability
- Height adjustable from 47" to 85" (119.4 cm to 216 cm)
- Base and pole made of stainless steel for durability
- Two ram hooks hold solution containers securely

ITEM# MDS80441
PKG. 1 EACH



HOMECARE

I.V. POLE STAND

Economical for light-duty home use

- Five-legged plastic base offers greater stability
- Base diameter of 23.5" (59.7 cm)
- Tool-free push-button design makes adjusting easy
- Height adjustable from 37.5" to 83.5" (95.3 cm to 212 cm)
- Two ram hooks hold items securely

ITEM# MDS129CS
PKG. 3/CS



OTHER POPULAR PRODUCTS

STAINLESS STEEL

ITEM# MPH08MS01C
PKG. 1 EACH



Mayo stand

- Tray measures 16" x 21" (41 cm x 53 cm)
- Foot operated for ease and convenience
- Height adjusts from 39.5" to 62" (100 cm to 157 cm)

ITEM# MPH08190
PKG. 1 EACH



Kickbucket

- 13 qt (12.3 liter)
- With wall guard rubber bumper and 4 roller blade wheels for smooth, wobble-free movement

ITEM# MPH08SS01C
PKG. 1 EACH



Stacking stool

- 12"L x 18"D x 6"H (30 cm x 46 cm x 15 cm)
- Stability bars for extra safe footing
- Connection points for secure assembly of multiple units

ITEM# MPH0807SSHB
PKG. 1 EACH



Single basin stand with H-brace

- 8.5 qt (8 liter) basin
- 15"L x 15.25"W x 33"H (38 cm x 39 cm x 84 cm)

GET ALL THE DETAILS

Contact your Medline sales representative to find out how we can help you with your equipment needs.



Medline Industries, Inc.
One Medline Place
Mundelein, IL 60060

Medline United States
1-800-MEDLINE (633-5463)

www.medline.com | info@medline.com

Medline Canada
1-800-396-6996

www.medline.ca | canada@medline.com

Medline México
01-800-831-0898

www.medlinemexico.com | mexico@medline.com